



## FOXX K1 Offroad compactbike

Approved up to 120 kg payload

### Basic equipment

- Front wheels 24" (507-30) tubless ready
- Rear wheels 27.5" (587-35) tubless ready
- Full suspension with 100 mm front travel and 140 mm rear travel
- A-shaped handles with integrated shift and brake controls
- Hydraulic disc brake system front and rear
- Backrest adjustable in angle and distance
- Bottom bracket position adjustable
- Mudguard for rear wheel
- Underride guard
- Battery lighting set

A9105100000

€ 16,975.00

☐ **Quotation date:** \_\_\_\_\_

**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

☐ **Order date:** \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

**Order no:** \_\_\_\_\_

### Medical suppliers / company address:

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**Contact:** \_\_\_\_\_

☐ **Different invoice address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

☐ **Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

### Customer name / commission:

**ID- / Com. no:** \_\_\_\_\_

**Surname, forename:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**D.O.B:** \_\_\_\_\_

**Paying organisation:** \_\_\_\_\_

### Body dimensions

**Body height:** \_\_\_\_\_ cm

**Upper body height:** \_\_\_\_\_ cm

**Weight:** \_\_\_\_\_ kg

**Indication:** \_\_\_\_\_

Information required to process the transaction.

## Drive system Bafang

The specified ranges depend on the route profile, user weight and speed.

### Bafang M560 motor

Maximum power for all off-road adventures

☐ A9105100010  
500W, 130Nm torque,  
assistance up to 25 km/h,  
starting assistance up to 15 km/h



### Battery Bafang 48V 696Wh

Range: approx. 50 km - 75 km

☐ A9105900100  
Battery mounted in the main frame under the  
seat plate



### Double battery Bafang 48V 696Wh

Range: approx. 100 km - 150 km

☐ A9105900105 **€ 855,00**  
Batteries mounted in the main frame under the  
seat plate.  
Interchangeable plug-in connection mounted  
on the seat plate for easy switching between  
batteries.



## Drive system Shimano

The specified ranges depend on the route profile, user weight and speed.

### Shimano EP6 Cargo motor

Drive system for off-road and on-road use

☐ A9105100020 **€ 325,00**  
250W, 85Nm torque,  
assistance up to 25 km/h,  
starting aid up to 6 km/h  
Pedelec system with road approval



### Battery Shimano 36V 630Wh

Range: approx. 45 km - 70 km

☐ A9105900110 **€ 610,40**  
Battery mounted in the main frame under the  
seat plate



### Double battery Shimano 36V 630Wh

Range: approx. 90 km - 140 km

☐ A9105900115 **€ 1,405.30**  
Batteries mounted in the main frame under the  
seat plate.  
Interchangeable plug-in connection mounted  
on the seat plate for easy switching between  
batteries.



## Colour attachments

If not specified, the legrest, rear mudguard and chain guard disc will be delivered in UDS black.

### Leg rest, rear mudguard and chain guard disc UDS black

☐ A9105100170-001

### Leg rest, rear mudguard and chain guard disc Pure orange RAL 2004 gloss

☐ A9105100170-002

### Leg rest, rear mudguard and chain guard disc Customized colour selection

☐ A9105100170-003 **€ 219,00**  
from colour palette one-colour design  
powder colours wheelchair

A9000100170- \_\_\_\_ \_

Colour: \_\_\_\_\_

## Derailleur gears

### Shimano Cues U8000

11-speed derailleur gear system

☐ A9105200010  
with Linkglide technology  
Cassette 11-50 T



### Sram GX-Eagle AXS

12-speed derailleur gear system

☐ A9105200011 **€ 767,50**  
electric radio control  
Cassette 11-50 T



## Gear hub

### Rohloff 14-speed gear hub with electric operation

☐ A9105200020 € 1,987.50

is supplied via the Bafang system.



## Chassis

### Standard chassis

☐ A9105400010  
Fox Float DPS Performance front and Manitou Mara Inline rear with three modes each (Open, Medium, Firm)  
Rebound and compression adjustable

### High-quality off-road chassis

☐ A9105400020 € 840,00  
Fox Float DPS Factory front and Manitou Mara Pro rear with three modes each (Open, Medium, Firm)  
Rebound and compression adjustable  
High- and low-speed compression adjustable

## Dimensions

### Crank type and grip width

#### Standard crank

- ☐ Grip width 40 cm  
A9117500100-040
- ☐ Grip width 45 cm  
A9117500100-045
- ☐ Grip width 50 cm  
A9117500100-050

#### Tetra crank (flared)

with extended crank arms for maximum upper body space

- ☐ Grip width 48 cm  
A9117500100-048
- ☐ Grip width 53 cm  
A9117500100-053
- ☐ Grip width 58 cm  
A9117500100-058

#### Standard Tretkurbel

- ☐ Length 155 mm  
A9110600100-155
- ☐ Length 175 mm  
A9110600100-175
- ☐ Length 195 mm  
A9110600100-195



#### Tetra Tretkurbel

- ☐ Length 155 mm  
A9100600105-155 € 132,10
- ☐ Length 175 mm  
A9100600105-175 € 132,10
- ☐ individual length  
140 mm to 195 mm  
A9100600100 € 251,20



### Adjustable bottom bracket support

adjustable in angle via clamped pivot point on the fork bridge and in length.

Please specify the bottom bracket post size that comes closest to your desired position for the presetting.

☒ S (180 mm)  
A9105500010-01

☐ M (240 mm)  
A9105500010-02

Please note: If not specified, the adjustable bottom bracket support will be delivered in S.



## Handles

### Palm rest on the handle

☐ A9110601200 € 99,80

upholstered, bolted, angle-adjustable.  
Not in combination with other handles.



### Tetra hand support, pivotable, adjustable grip distance

☐ A9110601101 € 164,50



### Tetra backrest swivel and adjustable grip width

☐ A9110601102 € 186,00

very stiff external support as round tube  
Ø 20 mm with foam rubber coating



## Upper body fixation

### Chest strap for upper body fixation

☐ size S (105 cm) € 231,80  
A9000800016

☐ size M (125 cm) € 231,80  
A9000800013

☐ size L (145 cm) € 231,80  
A9000800014

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.  
Colour: black



### Pelvi.Loc thigh fixation


☐ A9105800000 € 330,00



## Tyres


### Offroad tyres

☒ A9105105110  
Front: Smart Sam 24 x 2,1"  
Rear: Magic Marry 27,5 x 2,6" Ultra Soft



### Trekking tyres

☐ A9105105100 € 46,80  
Front: Marathon Plus 24 x 1,75"  
Rear: Johnny Watts 27,5 x 2,6"



### Tubless kit with tyre insert for rear wheel

☐ A9105105120 € 78,40

## Accessories

### Lighting for front and rear

☐ A9105900000 € 219,00

incl. reflectors for StVZO



Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 5 - 2025-04-08 / The prices are plus statutory value-added tax and shipping costs

☐ = Included in the basic price ☒ = Position selected

Order sheet FOXX K1 Offroad compactbike Rev. 00 UK

## LEGAL NOTICE - FOXX offroad handbike with SHIMANO or BAFANG drive system

PRO ACTIV offers the FOXX offroad handbike with a choice of two different drive systems. The utilisation options differ depending on the equipment and the applicable legal regulations.

### SHIMANO drive system (pedelec system according to EN 15194)

The FOXX with SHIMANO drive system fulfils the requirements of a pedelec in accordance with EU standard EN 15194. With this equipment, the vehicle is permitted for use on public roads and public paths in Germany without the need for an operating permit or an insurance licence plate. This also applies to use on cycle paths.

### BAFANG drive system (more powerful off-road system)

The FOXX with BAFANG M560 drive system has a motor output of 500 W and a torque of 130 Nm. PRO ACTIV expressly points out that the FOXX with BAFANG drive system may not be operated on public roads or on public paths, squares or cycle paths in Germany without an individual operating licence and without an insurance licence plate. In this case, operation is only permitted on private property and off-road terrain (e.g. on non-public trails or bike parks). Operation in public areas without a valid operating licence constitutes an infringement of the Compulsory Insurance Act (§ 1 PflVG) and can be prosecuted in accordance with § 6 Para. 1 PflVG.

Please also note: If you use the FOXX with BAFANG drive outside Germany, it is your responsibility to inform yourself in advance about the legal requirements and traffic regulations applicable in the respective country.

### Confirmation of knowledge

I confirm with my signature that I have read, understood and taken note of the above information. I undertake to use the FOXX with BAFANG drive system in Germany only with a valid operating licence on public roads and only on private property or off-road terrain without an operating licence. When using the vehicle in other countries, I will inform myself on my own responsibility about the legal requirements applicable there and comply with them.

---

**Customer name in block letters**

---

**Place, date**

---

**Signature (end customer)**