

Im Hofstätt 11  
D-72359 Dotternhausen  
Phone +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.com



## LIFT activ

### with mechanical seat height adjustment

- Manual seat height adjustment by 30 cm supported by individually adjusted lifting springs
- Lifting function triggered via ergonomic handle
- Rigid frame chassis firmly welded together
- Seat unit in rigid frame construction bolted to chassis
- Total weight starting at 14.5 kg

Maximum user weight limit 100 kg

A9065100000 € 7,880.00

☐ Quotation date: \_\_\_\_\_

Consultation was carried out by  
Staff at the medical supplies shop: \_\_\_\_\_

PRO ACTIV field representative: \_\_\_\_\_

☐ Order date: \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

### Medical suppliers/company address:

Cust. no.: \_\_\_\_\_

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Contact: \_\_\_\_\_

Shipping address: ☐ same as company address

☐ Other (see below)

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

### Customer name/commission:

Surname, \_\_\_\_\_

forename: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

D.O.B: \_\_\_\_\_

Paying \_\_\_\_\_

organisation: \_\_\_\_\_

Invoicing address: ☐ same as company address

☐ other (see below)

Cust. no.: \_\_\_\_\_

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Contact: \_\_\_\_\_

# Basic size

## Body size:

Body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg


Lower leg length: \_\_\_\_\_ cm  
(anatomic dimension: hollow of the knee to the foot sole):

Indication: \_\_\_\_\_

UL

This information is mandatory for further processing.

## Wheelchair dimensions:



Seat width <b>SB</b>	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46
----------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Seat depth <b>ST</b>	Short frame					Long frame				
	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48

Back height <b>RH</b>	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48				

Lower leg length <b>USL</b>	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	

Seat height back <b>SHh</b> (for 24" drive wheels)	<input type="checkbox"/> <u>41</u>	Possible on request:	<input type="checkbox"/> 42	<input type="checkbox"/> 43				
Seat height front <b>SHv</b> (for 24" drive wheels)	<input type="checkbox"/> <u>47</u>	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51	Possible on request:	<input type="checkbox"/> 45	<input type="checkbox"/> 46

Tipping point <b>KP</b> (front side rear tube to axis middle drive wheel)	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> <u>8</u>	<input type="checkbox"/> 9	<input type="checkbox"/> 10
--	----------------------------	----------------------------	-----------------------------------	----------------------------	-----------------------------

Please cross on the appropriate dimension (all dimensions are specified in cm).

Underlined values are standard settings.

The values SHh and SHv may deviate from the specified values depending on the wheel size and wheel camber!

Frame colour – one-colour design		
Information about current colour designs and example images at: <a href="http://www.proactiv-gmbh.de/farben.php">www.proactiv-gmbh.de/farben.php</a>		
frame colour: Selection from colour palette		
Frame colour is specified when ordering <input checked="" type="checkbox"/> A9000100170-000		
Powder coatings		
<b>White RAL 9016</b> <input checked="" type="checkbox"/> gloss A9000100170-023 <input type="checkbox"/> matt A9000100170-123 € 40.10	<b>Brilliant silver</b> <input checked="" type="checkbox"/> gloss A9000100170-011	<b>Anthracite RAL 7016</b> <input checked="" type="checkbox"/> gloss A9000100170-001 <input type="checkbox"/> matt A9000100170-101 € 40.10
<b>Black RAL 9005</b> <input checked="" type="checkbox"/> gloss A9000100170-003 <input type="checkbox"/> matt A9000100170-019 € 40.10	<b>Pure orange RAL 2004</b> <input checked="" type="checkbox"/> gloss A9000100170-022 <input type="checkbox"/> matt A9000100170-122 € 40.10	<b>Pink magenta RAL 4010</b> <input checked="" type="checkbox"/> gloss A9000100170-037 <input type="checkbox"/> matt A9000100170-137 € 40.10
<b>Pastel blue RAL 5024</b> <input checked="" type="checkbox"/> gloss A9000100170-036 <input type="checkbox"/> matt A9000100170-136 € 40.10	<b>Lemon green</b> <input checked="" type="checkbox"/> gloss A9000100170-005 <input type="checkbox"/> matt A9000100170-105 € 40.10	<b>Yellow RAL 1023</b> <input checked="" type="checkbox"/> gloss A9000100170-046 <input type="checkbox"/> matt A9000100170-146 € 40.10
<b>Red RAL 3001</b> <input checked="" type="checkbox"/> gloss A9000100170-047 <input type="checkbox"/> matt A9000100170-147 € 40.10	<b>Lagoon blue</b> <input checked="" type="checkbox"/> gloss A9000100170-052 <input type="checkbox"/> matt A9000100170-152 € 40.10	<b>Athletic grey</b> <input checked="" type="checkbox"/> gloss A9000100170-050 <input type="checkbox"/> matt A9000100170-150 € 40.10
<b>Muddy grey</b> <input checked="" type="checkbox"/> gloss A9000100170-051 <input type="checkbox"/> matt A9000100170-151 € 40.10		
“Dormant” powder coatings		
<b>Deep black, silver effect</b> <input checked="" type="checkbox"/> gloss A9000100170-007	<b>Sky blue</b> <input checked="" type="checkbox"/> gloss A9000100170-017	<b>Marine blue</b> <input checked="" type="checkbox"/> gloss A9000100170-045 <input type="checkbox"/> matt A9000100170-145 € 40.10
<b>Electric green</b> <input checked="" type="checkbox"/> gloss A9000100170-015 <input type="checkbox"/> matt A9000100170-115 € 40.10	<b>Violet</b> <input checked="" type="checkbox"/> gloss A9000100170-044	<b>Fire red</b> <input checked="" type="checkbox"/> gloss A9000100170-018 <input type="checkbox"/> matt A9000100170-118 € 40.10
Frame colour – “liquid coating” one-colour design		
€ 323.20		
<b>Storm (ice blue)</b> <input type="checkbox"/> gloss A9000100170-500 <input type="checkbox"/> matt A9000100170-501	<b>Coffee</b> <input type="checkbox"/> gloss A9000100170-520 <input type="checkbox"/> matt A9000100170-521	<b>Prairie (rose gold)</b> <input type="checkbox"/> gloss A9000100170-504 <input type="checkbox"/> matt A9000100170-505
<b>Pearl (mother of pearl-white)</b> <input type="checkbox"/> gloss A9000100170-524 <input type="checkbox"/> matt A9000100170-525	<b>Berry</b> <input type="checkbox"/> gloss A9000100170-526 <input type="checkbox"/> matt A9000100170-527	<b>Blue Ocean</b> <input type="checkbox"/> gloss A9000100170-522 <input type="checkbox"/> matt A9000100170-523
<b>Emerald (emerald green)</b> <input type="checkbox"/> gloss A9000100170-514 <input type="checkbox"/> matt A9000100170-515	<b>Mint</b> <input type="checkbox"/> gloss A9000100170-528 <input type="checkbox"/> matt A9000100170-529	<b>Chilli</b> <input type="checkbox"/> gloss A9000100170-516 <input type="checkbox"/> matt A9000100170-517

## Frame colour – two colour designs

Note: Frames with two colour designs in the gloss level of the design colour.

- ☐ **Two-colour design comprising one basic colour (on the back of the frame) and one design colour (on the front of the frame)** € 248.30

Please select from the following colours:

### Base colour (on the rear frame)

#### White RAL 9016

- ☐ A9000100170-200

#### Brilliant silver

- ☐ A9000100170-201

#### Anthracite RAL 7016

- ☐ A9000100170-202

#### Black RAL 9005

- ☐ A9000100170-203

#### Pastel blue RAL 5024

- ☐ A9000100170-206

#### Lemon green

- ☐ A9000100170-207

#### Yellow RAL 1023

- ☐ A9000100170-208

#### Lagoon blue

- ☐ A9000100170-212

#### Athletic grey

- ☐ A9000100170-210

#### Muddy grey

- ☐ A9000100170-211

### Design colour (on the front frame)

#### Black RAL 9005

- ☐ gloss A9000100170-259  
☐ matt A9000100170-260

#### Lemon green

- ☐ gloss A9000100170-250  
☐ matt A9000100170-251

#### Chilli

- ☐ gloss A9000100170-253  
☐ matt A9000100170-254

#### Pure orange RAL 2004

- ☐ gloss A9000100170-255  
☐ matt A9000100170-256

#### Coffee

- ☐ gloss A9000100170-278  
☐ matt A9000100170-279

#### Mint

- ☐ gloss A9000100170-286  
☐ matt A9000100170-287

#### Emerald (emerald green)

- ☐ gloss A9000100170-274  
☐ matt A9000100170-275

#### Storm (ice blue)

- ☐ gloss A9000100170-264  
☐ matt A9000100170-265

#### Prairie (rose gold)

- ☐ gloss A9000100170-266  
☐ matt A9000100170-267

#### Pearl (mother of pearl-white)

- ☐ gloss A9000100170-282  
☐ matt A9000100170-283

#### Berry

- ☐ gloss A9000100170-284  
☐ matt A9000100170-285

#### Blue ocean

- ☐ gloss A9000100170-280  
☐ matt A9000100170-281

## Colour - Attachment parts

Note: The attachment parts are always mounted in black.

### Drive wheel: Aluminium rim, hub and spokes

- ☐ Black A9000100171-003

### Caster fork and bearing block

- ☒ Black A9000100173-003

## Attachment parts colour design according to customer desires:

### Drive wheel: Aluminium rim and design inlet<sup>1</sup>

- ☐ Pure orange RAL 2004 A9000100171-022 € 193.10  
☐ White RAL 9016 A9000100171-023 € 193.10  
☐ Frame colour A9000100171-999 € 230.80  
☐ Selection of single colour coating (deviating from the frame colour) A9000100171-998 € 267.10  
Colour: \_\_\_\_\_  
A9000100170- ☐ ☐ ☐

### Caster wheel: 5" aluminium rim<sup>2</sup> (3-spoke design)

- ☐ Pure orange RAL 2004 A9000100174-022 € 44.40  
☐ White RAL 9016 A9000100174-023 € 44.40  
☐ Frame colour A9000100174-999 € 69.50  
☐ Selection of single colour coating (deviating from the frame colour) A9000100174-998 € 80.60  
Colour: \_\_\_\_\_  
A9000100170- ☐ ☐ ☐

## Special colour requests (frame or attachment parts)

### Frame colour individual procurement Selection from RAL colour palette

- ☐ A9000100170-100 € 383.10

### Degree of gloss

- ☐ gloss  
☐ matt

### Colour designation

\_\_\_\_\_

### RAL number

\_\_\_\_\_

**Note:** Attachment parts **not** possible in design colours

<sup>1</sup> Only in conjunction with Liteline hub (for drive wheel aluminium rim with 18 or 36 spokes, object).

<sup>2</sup> Only in conjunction with 5" caster wheel aluminium rim 3-spoke design easy running A9000501200 or good shock absorption A9000501210

## Frames

Note: Footrest area with V-shape  $x = SB - (0.25 \times USL + 6 \text{ cm})$


### Straight front end

☐ A9065100010

Knee angle 100°  
only with long frame

☐ A9065100011

Knee angle 105°




### Front end V-shape

☐ A9065100020

Knee angle 100°  
from SB 40 cm only possible in combination  
with long frame

☐ A9065100021

Knee angle 105°



## Push handles

### Back tube with integrated handles

☐ A9000401051-K € 120.60

straight

Also possible with lumbar curvature:

☐ 1 cm

☐ 2 cm

☐ 3 cm



### Aluminium push handles fixed in back tube

☐ A9000401050 € 120.60

from RH 23 cm  
from RH 42 cm with lumbar joint



### Push handles, horizontally screwed in back tube

☐ A9000401076 € 223.60

can be unscrewed  
from RH 37 cm with lumbar joint



### Safety push handles with continuous height adjustment

☐ A9000401060 € 200.10

Mount in the back tube  
from RH 30 cm with adaptable back  
Starting at a back height of 35 cm for back  
tube with lordosis support  
from RH 42 cm with lumbar joint  
from RH 43 cm with back shell



### Safety push handles back-positioned length 45 cm

☐ A9000401075 € 344.80

continuous height adjustment,  
Safety push handle holder mounted on the  
back cross bar  
Up to seat width of 42 cm or body weight of 95  
kg



## Backrest

### Folding backrest

☒ A9065400011

7-way angle adjustment,  
With integrated plastic coated aluminium side sections

Note: The backrest can only be folded down completely from SB 41 cm. Up to SB 40 cm, we therefore recommend the additional option "back tube with lumbar joint".



### Back tube with lumbar curvature

☐ A9000400005-K € 105.80

☐ 1 cm

☐ 2 cm

☐ 3 cm



### Back tube with lumbar joint

☐ A9065401020-K € 391.00

Backrest can be additionally folded down at the level of the lumbar support from RH 35 cm

☐ 1 cm

☐ 2 cm

☐ 3 cm



## Back equipment

### Adaptable back with back padding

☒ A9000802010-K

Cover: Black, breathable, open surface  
Equipment with back cross bar



### Raised back + 3 - 5 cm

☐ A9000802100 € 144.40

for the nominal back height dimension  
Only in combination with height adjustable push handles



### Equipment without back upholstery and back padding or back system

☐ A9000802011 € - 53.00

Equipment with back cross bar


Other back systems, such as ergonomic back shell etc. see separate order sheet: "PRO ACTIV seat and back systems"



## Seat upholstery


### Body Contour

☒ A9000801025-K



### Open belt system

☐ A9000801010-K € 71.10





### Equipment without seat upholstery

☒ A9000801050



## Seat cushion



### Seat cushion in rectangular contour Inlet: Standard

☐ A9000805180-K\* € 106.60

### Seat cushion with radius contour Inlet: Standard

☐ A9000800300-K\* € 148.00

### Seat cushion height:

☐ 3 cm ☐ 5 cm  
☐ 4 cm ☐ 6 cm

### Seat cushion in rectangular contour Inlet: Anti-Decubitus

☐ up to 30x30 cm  
A9000805204-K\* € 276.20

☐ >30x30 up to 50x50  
cm  
A9000805200-K\* € 401.90

### Seat cushion with radius contour Inlet: Anti-Decubitus

☐ up to 30x30 cm  
A9000800311-K\* € 326.50

☐ >30x30 up to 50x50  
cm  
A9000800310-K\* € 443.70

### Seat cushion height:

☐ 4 cm  
☐ 6 cm

### Cover:

☐ Black, breathable, open surface

☐ Black, breathable, closed surface

If not specified, the cover fabric "black, breathable, open surface" is delivered.

☐ **Seat cushion is existing**

Brand/type \_\_\_\_\_

Cushion height \_\_\_\_\_ cm

You can also find the seat cushions on our separate "PRO ACTIV Seat and Back Systems" order sheet

## Lift mechanism

### Mechanical seat height adjustment by 30 cm

via gas springs up to approx. 100 kg

☐ Release on the right  
A9065300401-K\*

☐ Release on the left  
A9065300402-K\* € 102.40

Note: The release mechanism is only possible  
on the right **OR** on the left

## Clothing guard

### Aluminium clothing guard

rigidly attached to the armrest  
cover following the wheel arch



in anthracite matt

☐ A9065403210-K\*

in frame colour (basic colour with two-colour  
design)

☐ A9065403211-K\* € 74.10

in design colour (only possible with frame  
colour in two-colour design)

☐ A9065403212-K\* € 111.30

## Wheel camber

\*Total width with tight-mounted handrim.

### 1° drive wheel camber

☐ A9065202010-K  
Total width = approx. SB + 21 cm\*  
from SB 38 cm

### 2.5° drive wheel camber

☐ A9065202020-K  
Total width = approx. SB + 23.5 cm\*  
from SB 36 cm

### 4° drive wheel camber

☐ A9065202030-K  
Total width = approx. SB + 26 cm\*



## E-adaptations

### Preparation for "e-motion" mounting

☐ A9065201200\* € 465.00

Equipment: e-motion torque support,  
1/2" axle and socket length 50.8 mm (2")  
4° camber only up to 60 kg

Wheelchair delivery to

☐ Alber GmbH  
(Please also enter on p.1 under "Shipping  
address, other")

☐ Specialist dealer

## Brake system

### Brake system push to lock Standard brake handle

☐ A9000701010-001



### Brake system push to lock Long brake lever

☐ A9000701010-002



### Integrated parking brake

☐ A9065702020 € 229.80



### Brake lever at 20°

☐ A9000700002  
Only in conjunction with knee lever brake  
system, standard brake lever, long or low.



## Footrest

### Continuous footrest

☐ A9016601010-K\*  
with tubular resting bracket  
height adjustable +/- 3 cm  
For equipment of the footrest with "Rubber  
coated footplate support", please select under  
"Footrest accessories"



### Footbar folds up to the rear with spring locking mechanism

☐ A9016602005-K\* € 265.70

with continuous tubular resting bracket  
from USL 40 cm  
only possible with V-shape and frame long  
frame  
For equipment of the footrest with "Rubber  
coated footplate support", please select under  
"Footrest accessories"



### Footrest folding back to the rear

☐ A9016602010-K\* € 284.70

with continuous tubular resting bracket and  
20 mm operating lever  
from USL 42 cm  
with V-shape from SB 38 cm  
For equipment of the footrest with "Rubber  
coated footplate support", please select under  
"Footrest accessories"



### Divided footrest

☐ A9016603200-K € 311.40

with rubber coated footplate,  
height adjustable +/- 3cm, angle adjustment,  
can be pivoted outwards and folded up to the  
side  
including frame reinforcement  
with V-shape from SB 40 cm



### Removable frame stem, foldable with divided footrest

☐ A9065630000-K € 782.00

height adjustable, with rubber plate  
Knee angle 105°  
with straight front end only  
with short frame only with integral brake





## Footrest accessories

### Rubber coated footplate

☐ A9000610000-K € 116.40

bolted, angle-adjustable  
Depth: 14 cm



### Footplate side stop

☐ A9000610100 € 55.20

adjustable  
only in combination with rubber coated  
footplate



## Adaptations

### Anti-tipping support, pivotable upwards

☒ A9065300100  
removable using Quick-pin



## 4" caster wheel

### 4" caster wheel "Standard" solid rubber black

☐ A9000501011-004  
(82A), puncture-proof



### 4" Skate wheel (100 mm)

☐ A9000501120 € 65.40  
puncture-proof



### 4" safety light-up wheel (100 mm)

☐ A9000501220 € 87.60  
puncture-proof



### 4" Frog Legs caster wheel (100 mm)

☐ A9000501230 € 196.50  
with replaceable rubber 6 mm, puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



## 5" caster wheel

### 5" caster wheel "Standard" PU black

☐ A9000501036-005  
(82A), puncture-proof



### 5" Caster wheel "Standard" PU grey









☐ A9000501031-005  
(82A), puncture-proof



### 5" Caster wheel "2 components"

☐ A9000501030 € 52.00  
puncture-proof



<b>5" safety light-up wheel (125 mm)</b> <input type="checkbox"/> A9000501221 € 87.60 puncture-proof 	<b>5" Aluminium rim 3-spoke design "easy running"</b> <input type="checkbox"/> Silver A9000501200 € 160.50 <input type="checkbox"/> black A9000501201 € 160.50 With exchangeable rubber ring 3 mm hard, puncture-proof Incl. hair guard spacer (3 mm) 	<b>5" aluminium rim 3-spoke design "Good shock absorption"</b> <input type="checkbox"/> Silver A9000501210 € 166.30 <input type="checkbox"/> black A9000501211 € 166.30 with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm) 
<b>5" Frog Legs caster wheel (125 mm)</b> <input type="checkbox"/> A9000501232 € 196.50 with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm) 		
<b>6" castor wheels, on request after testing, only with V-shape front end</b>		
<b>Caster fork</b>		
Caster fork colours, see "Component colours"		
<b>Caster fork Standard</b> <input type="checkbox"/> A9000500100-K001 with screwed axle <input type="checkbox"/> A9000500100-K002 € 154.80 with quick-release axle 	<b>Heavy duty caster fork for high loads</b> Extra weight approx. 100 - 150 g <input type="checkbox"/> A9000510111-K001 € 49.90 with screwed axle <input type="checkbox"/> A9000510111-K002 € 202.40 with quick-release axle 	
<b>Caster fork spare wheel "Standard"</b> <input type="checkbox"/> A9000500300-K € 274.40 Standard caster forks with quick-release axles 	<b>Caster fork spare wheel "Heavy Duty"</b> <input type="checkbox"/> A9000500301-K € 312.50 Heavy duty caster forks with quick-release axles 	The spare caster fork does not have a caster wheel. In case a caster wheel is requested, please put a cross on a second selection in the category "Caster wheels".

## Drive wheels

Information on further rim types, versions and availabilities under: [www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)

### Drive wheel, 18 spoke aluminium rim

☐ 24" (ETRTO 540)  
A9000203122

with Liteline hub, spoked radially  
Not in combination with other puncture-proof  
tyres  
only up to a body weight of 70 kg:

Sturz \ Gewicht	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	✓
61 - 70 Kg	✓	✓			



### Drive wheel, 36 spoke aluminium rim

☐ 24" (ETRTO 540)  
A9000203000

Up to 120 kg body weight with Liteline hub  
(left image)  
Over 120 kg body weight with aluminium high  
flange (right image)  
radially spoked



### Drive wheel double rim 36 spokes

☐ 24" (ETRTO 540) € 145.90  
A9000203100

☐ 25" (ETRTO 559) € 145.90  
A9000203200

☐ 26" (ETRTO 590) € 145.90  
A9000203300

with aluminium front hub,  
radially spoked  
Use with high loads



### Spinergy running wheel "Spox Sport" with 24 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540) € 725.00  
A9000203550

☐ 25" (ETRTO 559) € 725.00  
A9000201559

☐ 26" (ETRTO 590) € 725.00  
A9000203555



### Spinergy "Everyday Light Extreme" with 12 PBO spokes; rim and hub in black, spokes white

☐ 24" (ETRTO 540) € 885.00  
A9000203554




### Spinergy "Everyday Light Extreme" with 12 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540) € 885.00  
A9000203552



## Handrims

x = dimension handrim attachment



Handrim attachment

General:

☐ x = 11 mm

☐ x = 19 mm

Handrim attachment

Curve handrim:

☐ x = 10 mm

☐ x = 18 mm

Unless otherwise specified, the handrim is mounted wide (19 mm / 18 mm)

### Stainless steel Ø 20 mm

☐ 24" handrim for  
24" wheel and 25"  
wheel  
A9000205010-004

☐ 26" handrim for  
26" wheel  
A9000205010-006



### Aluminium Ø 19 mm silver anodised

☐ 24" handrim for  
24" wheel and 25"  
wheel  
A9000205019-004



### Aluminium Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

☐ = Gloss surface

☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

☐ 24" handrim for  
24" wheel and 25"  
wheel  
A9000205020-004 **€ 63.10**

☐ 26" handrim for  
26" wheel  
A9000205020-006 **€ 63.10**



### Aluminium Ø 19 mm polished and hard-anodised

☐ 24" handrim for  
24" wheel and 25"  
wheel **€ 152.50**  
A9000205030-004

☐ 26" handrim for  
26" wheel **€ 152.50**  
A9000205030-006



### Titanium Ø 19 mm

☐ Handrim for  
24" wheel and 25"  
wheel **€ 310.30**  
A9000205040-004

☐ 26" handrim for  
26" wheel **€ 310.30**  
A9000205040-006









### Silicone handrim cover

☐ for 24" handrim  
A9000205100-004 **€ 215.20**

☐ for 26" handrim  
A9000205100-006 **€ 215.20**

Surface without structure,  
compatible with handrim Ø 19 mm or  
Ø 20 mm.



<p><b>Ultra-Grip “very grippy”</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205070-004 <b>€ 372.10</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel A9000205070-006 <b>€ 372.10</b></p> 	<p><b>Ergo-Grip</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205080-004 <b>€ 377.90</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel A9000205080-006 <b>€ 377.90</b> Attachment fitting to the rim</p> 	<p><b>Curve black anodised, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 24 mm A9000205130-004 <b>€ 370.60</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L A9000205140-004 <b>€ 370.60</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel H 31 mm, size L A9000205140-006 <b>€ 370.60</b></p> 
<p><b>Curve silver anodised ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L A9000205160-004 <b>€ 375.50</b></p> 	<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = Frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 24 mm A9000205135-004 <b>€ 432.40</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L A9000205145-004 <b>€ 432.40</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel H 31 mm, size L A9000205145-006 <b>€ 432.40</b></p> 	<p><b>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm size L A9000205150-004 <b>€ 481.40</b></p> 

<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205170-004 <b>€ 492.90</b></p> <p><input type="checkbox"/> 25" handrim for 25" wheel A9000205170-005 <b>€ 492.90</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel A9000205170-006 <b>€ 492.90</b></p> 	<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205175-004 <b>€ 477.30</b></p> <p><input type="checkbox"/> 25" handrim for 25" wheel A9000205175-005 <b>€ 477.30</b></p> <p><input type="checkbox"/> 26" handrim for 26" wheel A9000205175-006 <b>€ 477.30</b> with a robust and replaceable silicone lip, fully integrated into the profile</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205185-004 <b>€ 488.80</b> with a robust and replaceable silicone lip, fully integrated into the profile</p> 
<p><b>Gekko S</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205195-004 <b>€ 477.30</b> black anodised, with complete integration of the integrated, robust and replaceable silicone lip into the profile Profile height 24 mm</p> 	<p><b>Equipment without handrims</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>	
<p><b>Spoke guard</b></p> <p><b>Finger- and spoke guard "transparent"</b></p> <p><input type="checkbox"/> for 24" wheel A9000206300-004 <b>€ 191.50</b></p> <p><input type="checkbox"/> for 26" wheel A9000206300-006 <b>€ 191.50</b></p> 		

## Tyres

### Schwalbe "RightRun" High-pressure tyre

- ☐ 24x1" (540-25)  
A9000204001-004  
Colour: black  
7.5 bar



### Schwalbe "RightRun" High-pressure tyre

- ☐ 24x1" (540-25)  
A9000204000-004  
☐ 25x1" (559-25)  
A9000204000-005  
☐ 26x1" (590-25)  
A9000204000-006  
Colour: Grey, black  
7.5 bar



### Schwalbe "DownTown" Profiled tyres

- ☐ 24x1" (540-25)  
A9000204010-004  
7.5 bar



### Schwalbe "Marathon Plus Evolution" with "smart guard" puncture protection

- ☐ 24x1" (540-25) € 45.50  
A9000204030-014  
☐ 25x1" (559-25) € 45.50  
A9000204030-015  
☐ 26x1" (590-25) € 45.50  
A9000204030-016  
Colour: black  
10 bar



### Schwalbe "Marathon Plus" with "smart guard" puncture protection

- ☐ 24x1" (540-25) € 45.50  
A9000204030-004  
Colour: Grey, black wall  
10 bar



### Schwalbe "One"

- ☐ 24x1" (540-25) € 86.00  
A9000204123-004  
Colour: black  
9 bar



### Puncture-proof tyre "solid rubber tyres"

- ☐ 24x1" (540-25) € 113.40  
A9000204201-004  
Colour: black  
Not in combination with Spinergy drive wheels



## Accessories

### Through axle with assistance, Large push button

- ☐ A9000200096 € 97.80  
Interlock dimension: 122 / 110 mm  
Axle diameter: 12 mm  
Not in conjunction with lockable drive wheel  
axle clip "Tetra Clip",  
Not in combination with one-hand drive



### Lockable drive wheel axle clip "Tetra Clip"

- ☐ A9000206400 € 153.90  
Not in combination with Spinergy drive  
wheels, one-hand drive, drum brake or  
Liteline hub



### Hand air compressor with pressure gauge

- ☐ A9000901100 € 98.30  
up to 14 bar





<b>Length-adjustable safety belt with buckle</b> <input type="checkbox"/> A9000902100      € 74.60 	<b>Adjustable safety belt with hook and fastener</b> <input type="checkbox"/> A9000902110      € 59.60 	<b>Crutch holder with nylon cup</b> <input type="checkbox"/> A9000902310      € 127.00 flexible 
<b>SIGMA lighting set</b> <input type="checkbox"/> A9000903270      € 99.80 Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations 	<b>Transport bag for drive wheels</b> <input type="checkbox"/> for 24" wheel A9000903530-004      € 71.90 Ø 76 cm x 2 cm <input type="checkbox"/> for 26" wheel A9000903530-006      € 71.90 Ø 83 cm x 2 cm 	<b>Accessory bag</b> <input type="checkbox"/> A9000903600      € 70.30 17 cm x 13cm x 5 cm 
<b>Further accessories, e.g., adaptive bikes, steering attachments, etc. see separate order sheet or <a href="http://www.proactiv-gmbh.com">www.proactiv-gmbh.com</a></b>		

Notes and remarks: