



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Different invoicing address (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____
Indication: _____

FREAK compact bike

approved up to 100 kg payload

- Underframe chassis with reinforcement strut
- chassis length 102 cm, SB 35 cm
- Drive unit can be uncoupled via excenter lock
- black drive unit with 20" drive wheel
- A-shaped handles with integrated gearshift controls
- Shimano V rim brake
- Backrest with continuous angle adjustment
- wheel camber 6°, rear wheel size 20"
- safety flag
- Seat upholstery Body Contour

A9113000013 € 4,640.00

Drive wheel on Freak with 26"

A9113000014 € 253.60

V-brake replaced by mechanical disc brake.
Not compatible with Nexus 8-speed hub.

Body size

Body height: _____ cm
Weight: _____ kg

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

white RAL 9016

- gloss A9000100170-023
 matt A9000100170-123 € 43.40

brilliant silver

- gloss A9000100170-011

anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101 € 43.40

black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019 € 43.40

deep black effect silver

- gloss A9000100170-007

pure orange RAL 2004

- gloss A9000100170-022

fire red

- gloss A9000100170-018
 matt A9000100170-118 € 43.40

pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137 € 43.40

sky blue

- gloss A9000100170-017

pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136 € 43.40

navy blue

- gloss A9000100170-045
 matt A9000100170-145 € 43.40

electric green

- gloss A9000100170-015
 matt A9000100170-115 € 43.40

lemon green

- gloss A9000100170-005
 matt A9000100170-105 € 43.40

violet

- gloss A9000100170-044

yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146 € 43.40

red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147 € 43.40

lagoon blue

- gloss A9000100170-052
 matt A9000100170-152 € 43.40

athletic grey RAL 7046

- gloss A9000100170-050
 matt A9000100170-150 € 43.40

muddy grey RAL 7033

- gloss A9000100170-051
 matt A9000100170-151 € 43.40

peach RAL 3022

- gloss A9000100170-038
 matt A9000100170-138 € 43.40

Neodrives drive system

e-assistant version with neodrives

- A9113200000 € 3,250.10

- Neodrives Z10 motor
 - torque 12/40 Nm, power 250 W
 - colour display sMMi Z15 Connect
 - push assist for riding without pedal assistance up to 6 km/h
 - Neodrives RR G2 battery with 14.25 Ah / 513 Wh
 - range approx. 40 - 70 km
- Not in combination with hub gears!

>> Expected to be available from August 2026 <<



Hub gear

Not compatible with e-assistant version.

Shimano „Nexus“ 8-speed hub gear

- A9113000015
 „Grip Shift“ twist shift
 Gear ratio: 32 teeth sprocket,
 cassette 20 teeth
 with back-pedal brake
 Shimano V-Brake brake handle left

Shimano „Alfine“ 11-speed hub gear

- A9113200026 € 833.70
 „Rapid Fire“ gear lever
 Gear ratio: 32 teeth sprocket,
 cassette 20 teeth
 with hydraulic disc brake (Magura MT2) brake
 handle right and Shimano V-Brake brake
 handle left



Derailleur gears

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

A9113200034 € 517.50

mechanical derailleur gear with control on the handle
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT4) brake lever on the right
and Shimano V-brake brake lever on the left
Chainring at the top with chain guard discs

Selection of the chainring:

- 32 teeth
 38 teeth

Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

A9113200035 € 1,800.20

electronic derailleur gear with operation on the handle.
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT4) brake lever on the right
and Shimano V-brake brake lever on the left
Chainring at the top with chain guard discs

Selection of the chainring:

- 32 teeth
 38 teeth



Brake system

Powerful braking system with brake disc Ø 203 mm

A9110300105 € 83.10

recommended from 80 kg body weight.



Dimensions

Crank type and grip width

Bottom bracket shaft with standard crank

- grip width 40 cm
A9110500100-040
 grip width 45 cm
A9110500100-045

Bottom bracket shaft with tetra crank (flared)

- with extended crank arms for maximum
upper body space
 grip width 48 cm (176 mm)
A9110500100-048

Standard pedal crank

- length 155 mm
A9110600100-155
 length 175 mm
A9110600100-175



Tetra pedal crank


- length 155 mm € 136.10
A9100600105-155
 length 175 mm € 136.10
A9100600105-175



Fixed bottom bracket support

(Limited availability – only as long as stocks last. Alternative: Adjustable bottom bracket support)

- S (X/Y 67/153 mm)
A9100501100-003



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

- S (180 mm) A9100500001-002 € 232.10
 M (240 mm) A9100500001-003 € 232.10

If not specified, the adjustable bottom bracket support will be delivered in M.



Legrest

Divided legrest

- A9100100008
adaptation on the fork
for short legs
Only possible up to body size 165 cm



Foot stop for legrest

- A9100100012 € 91.20



Back equipment

Backrest with filling plate and back cushion

(Limited availability – only as long as stocks last.)

- A9113800002
Cover: Black, breathable, open surface

Backrest upholstery Body Contour


- A9100800010 € 92.40
with back cushion 3 cm
Cover: Black, breathable, open surface



Seat cushion

Seat cushion „standard“

- A9113000190 € 115.00
for chassis 102 cm
SB: 35 cm
height: 4 cm
Cover: Black, breathable, open surface



Tyres

20x1.1" Schwalbe „Durano“

A9113105001
(ETRTO 406-28)
road bike tyres



20x1.75" Schwalbe „Marathon Plus“

A9113105120 € 71.10
(ETRTO 406-47)
puncture-proof
with Smart Guard puncture protection inserts



Accessories

Plug-in mudguard

A9100502100 € 58.10



Rear guard clamped on the axle tube

A9100800100 € 218.10
with 40 mm mounting clamp



Neck rest

A9100800150 € 261.00
height-adjustable, angle-adjustable



Pelvic fixation belt

A9100800200 € 121.70



SIGMA BC 8.0 ATS bike computer

A9110901000 € 70.10
wireless, with mount
Not in combination with e-assistant version



SIGMA BC 14.0 WL CAD bike computer

A9110900018 € 171.50
with pedal frequency set, wireless
Not in combination with e-assistant version



Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

A9110903271 pcs. á € 23.20

Colour: Natural anodised

_____ piece



SIGMA lighting set

A9120900019 € 123.90

consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable.
Lighting set compliant with road traffic
regulations



Options for users with reduced functions

Palm rest on the handle

A9110601200 € 102.80
upholstered, bolted, angle-adjustable.



Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 169.40



Tetra backrest swivel and adjustable grip width

A9110601102 € 191.60

very stiff external support as round tube
Ø 20 mm with foam rubber coating
Not in combination with Tetra hand support
pivotable



Tetra shell handles A9110601100 € 172.40**Services****Consultation, clarification of the installation situation and measurement through PRO ACTIV** A9990000000 € 354.50**Handover and instruction of the handbike / traction device by PRO ACTIV** A9990000011 € 354.50**Notes and remarks:****Location, Date**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 6 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs

= included in the basic price = Position selected

Order sheet FREAK compact bike Rev. 07 UK