



WHEEL-E

up to 150 kg trailer load and 10 kg payload,
equipped with:

- Docking plate for quick connection to the „Quick-Snap adapter”
- Display for driving speed, assistance level and charge level display
- Driving speed up to 6 km/h with back gear, adjustable via drive signal transmitter on the right
- with hydraulic disc brake system and additional mechanical disc brake with parking function
- Parking stand with maneuvering rollers
- Handlebar, attached to the angle-adjustable stem
- Sigma LED battery lighting set with front and rear light
- Bell
- Drive wheel 20" (ETRTO 406), rim black
- Drive unit black

„The powerful, high-performance drive for demanding conditions”

Climbing power up to 15%

A9130100031-020

€ 6,530.00

Standard motor

55 Nm torque / rated power 1000 W

☐ included in the basic price

Driving comfort control with a more powerful HT motor

70 Nm torque / rated power 1000 W

☐ A9130100050

€ 496.30

incl. downhill brake in the speed control function with active motor brake. One-hand operation possible. Accelerate and brake with the thumb throttle.

Körpermaße

Body height: _____ cm

Upper body height: _____ cm

Weight: _____ kg

This information are required for the processing.

☐ Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

☐ Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

☐ Different invoice address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

☐ Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying organisation: _____

Indication: _____

Brake system

**parking disk brake,
operating handle mounted on the
handlebar support tube**

☐ A9130300106 € 156.50

designed as a tension spreader.
Not in conjunction with double disk brake
system

**Double disc brake system „Magura“
brake disc Ø 203 mm**

☐ A9130300105 € 197.10

recommended from 80 kg trailer load
and steep gradients, brake locking function via
locking bracket



Batteries and accessories

Battery 48V 10,5Ah 504Wh

☐ A9130900016
Range: approx. 25-45 km
Battery mounted on the right, balancing weight
mounted on the left



Double Battery 48V 10,5Ah 504Wh

☐ A9130900017 € 1,513.20

Range: approx. 50-90 km
with plug-in connector positioned on the
handlebar via magnetic plug
incl. 2nd charger



4.3" LCD color display with USB port

☐ A9130200030 € 225.70

mounted on the right side of the handlebar
clamp (replaces the 2.6" standard display)

Activation max. speed 15 km/h

☐ A9130900101 € 140.80

Use only on private property, outside the scope
of application of road traffic regulations

Activation max. speed 24 km/h

☐ A9130900106 € 140.80

Use only on private property, outside the scope
of application of road traffic regulations

LEGAL WARNING for the WHEEL-E towing device.

Equipping the WHEEL-E with a console for a maximum speed of 6 km/h.
The factory settings of the WHEEL-E PRO ACTIV are designed for an electric drive system that
supports a maximum speed of 6 km/h.
Consoles for a top speed of 15 km/h or 24 km/h are optionally available. PRO ACTIV points out
that the specific laws of your country must be observed for the use of the 15 km/h and 24 km/h
consoles. When ordering, please confirm with your signature that you have taken note of this fact.

Tyres

20x2.35" Schwalbe „Smart Sam“

☐ A9117105101
(ETRTO 406-60)
Profile tyres
Puncture protection level 3.5



20x2.40" Schwalbe „Super Moto-X“

☐ A9130105132 € 27.90
(ETRTO 406-62)
Puncture protection level 6



20x4.0" BFT tyres


☐ A9130105130 € 42.70
(ETRTO 406-100)
Profile tyres, Fat Tire



Handlebar

Steering column

- ☐ A9130600001
length of the steering column 220 mm
- ☐ A9130600002
length of the steering column 280 mm
- ☐ A9130600003
length of the steering column 340 mm
angle adjustable via clamped pivot point on the fork bridge



Accessories

SKS mudguard

- ☐ A9130500003
removable



Mudguard MX Colour: orange

- ☐ A9130500002
fixed mounted



Mudguard MX Colour: Black

- ☐ A9130500004
fixed mounted



Aluminum luggage rack „cargo systems”

- ☐ A9110903290 **€ 183.30**
Color: Black
Max. load weight 10 kg on the loading platform,
Max. 15 kg in conjunction with luggage bags at
lower loading height.



Bottle holder incl. drinking bottle

- ☐ A9130900001 **€ 23.20**
mounted on the left side of the battery
holder rail
Not in conjunction with double battery



LED front headlight

- ☐ A9130900004 **€ 118.30**
Power supply for the LED front headlight via
the drive system battery. Rear light as plug-in
light with battery supply. Lighting set omitted.

Rear-view mirror left

- ☐ A9130900002 **€ 90.60**

Rear-view mirror right

- ☐ A9130900003 **€ 90.60**

Options for users with reduced functions

Twist grip Triride

☐ A9130200019 € 45.30

instead of the thumb throttle handle the twist grip is installed.
Mounted on the right side of the handlebar.



Gas assistance for WHEEL-E

☐ A9130200017 € 27.20

tool-free clamping on the twist grip of the WHEEL-E



Gas operation Tetra

☐ A9130400100 € 334.40

right side



Brake operation Tetra

☐ A9130400101 € 364.90

brake operating on the left,
with brake operation
via rotary function with hand rest



Chest strap for upper body fixation

☐ size S (105 cm) A9000800016 € 231.80

☐ size M (125 cm) A9000800013 € 231.80

☐ size L (145 cm) A9000800014 € 231.80

☐ size XL (165 cm) A9000800017 € 231.80

consisting of elastic webbing material in
combination with neoprene with Velcro and
operating tab.
Colour: black



Quick release clip lever extension for rigid frame adapter

☐ A9110700004 € 100.90

to make the opening and closing of the quick
release clips easier.
only in conjunction with rigid frame adapter



Excenter lock handle extension for Quick Snap connection


Simplifies the opening and closing of the
Quick Snap adapter. Tetra version

☐ A9110700005 € 71.70

straight version
only in conjunction with rigid frame adapter

☐ A9110700006 € 71.70

curved version
only in conjunction with front adapter



Adaptation help for centering pin

☐ A9110700002 € 53.40

to simplify the insertion of the central adapter.



Parking stand with operating handle for high comfort of use

☐ A9110900004 € 57.00

with longer operating handle

Frame adapter

SPEEDY rigid frame adapter with "Quick-Snap-Adaption"

with quick-release clamp and centering pin holder

Colour: black

☐ A9110700141 € 502.40

Heavy duty version

☐ A9110700143 € 739.90

up to a max. body weight of 200 kg



TRAVELER folding frame adapter with "Quick-Snap-Adaption"

☐ A9110701300 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism

Colour: black



BUDDY rigid frame adapter with "Quick-Snap-Adaption"

☐ A9110701400 € 502.40

with quick-release clamp and centering pin holder

Colour: black



Wheelchair front adapter with universal clamp and "Quick-Snap-Adaption"

☐ A9110701500 € 572.30

adapter attachment to the wheelchair stem

Please select below:

☐ A9110701525

Tube diameter 23-27 mm

Mounting situation for PA products:

Straight front end from SB 31 cm to 54 cm

V-shaped front end from SB 30 cm to 50 cm

PLEASE NOTE:

If the exact tube diameter is not specified, the standard adapter for a tube diameter of 23-27 mm will be delivered.

☐ Clamp sleeve [70mm] included in standard


For close distance between footrest and drive wheel

☐ Clamp sleeve [110mm] € 44.50
A9199000117

For medium distance between footrest and drive wheel

☐ Clamp sleeve [170mm] € 53.90
A9199000154

For a wide distance between footrest and drive wheel



Clamp lever with 3-way star handle for actuating the frame clamping

☐ A9199701507 € 103.00

mounted on the right and left for easy terminal operation with limited manual function

Only in conjunction with front adapter



Accessories depending on wheelchair

Extended bike plate

- ☐ For SPEEDY A1 / A2
A9110100100-002 € 359.20
- ☐ For BUDDY
A9110100100-003 € 359.20
- ☐ For TRAVELER
A9110100100-001 € 359.20
- 4° camber integrated
Track independently adjustable



Wheelbase extension via „locking mechanism with pull operation”

- ☐ For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you
A9110102001 € 548.30
- removable, with continuous axle tube
4° / 6° camber integrated
Track independently adjustable



Wheelbase extension angle-adjustable and removable

- ☐ For SPEEDY F2 / F4
A9015100005 € 711.20
- with continuous axle tube,
height adjustment via toothed joint
adjustable in 6° angle increments
Wheel position displacement approx. 14 cm
4° / 6° camber integrated, toe adjustable
Mounting option must be checked by PRO
ACTIV.
Only possible in conjunction with base support
plate.
up to 120 kg trailer load.



Front brake

- ☐ A9125500120 € 148.10
- as hill start assist with anti-roll back device,
compatible for installation in combination with
the wheelchair knee lever brake.
Required for WHEEL-E with 15 km/h individual
operating permit.
- For combination with the wheelbase extension,
the article A9000700000 / A9000700173
holding rails adjustment range 295 mm is also
required.



Services

Delivery and instruction of the handbike / traction device by the specialist dealer incl. assembly of the frame adapter

on a PRO ACTIV wheelchair

☐ A9990000008

prod. no: _____

on a third-party product

☐ A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks: