

Im Hofstätt 11
D-72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



BUDDY classic **Active easy running wheelchair**

Maximum user weight limit 100 kg

A9041100000* € 3,690.00

☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

PRO ACTIV field representative: _____

☐ **Order date:** _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address: ☐ **same as company address**

☐ **other** (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: ☐ **same as company address**


☐ **Other** (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name/commission:

Surname, _____
forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____

Body size:



Body height: _____ cm

Weight: _____ kg



Lower leg length: _____ cm

(anatomic dimension: hollow of the knee to the foot sole):

Indication: _____

This information is mandatory for further processing.

Abducted frame shape, 2.5 cm/side



Unterschenkelänge **USL**
(Rollstuhlmaß unter Berücksichtigung
der Sitzkissenhöhe)

Seat width SB	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
	The following dimensions are optionally available subject to surcharge:				<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46

Size M (Frame length 34 cm (to the bend) and frame height 21 cm)											
<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31		

Seat depth **ST**

Size L (Frame length 38 cm (to the bend) and frame height 21 cm)											
<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35		

Size XL (Frame length 43 cm (to the bend) and frame height 21 cm)											
<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40		

Back height RH	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				
	Optional in size XL:			<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48				

Lower leg length **USL**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44								
Optional in sizes L + XL:			<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47					

Seat height back SHh	for 20" drive wheel						
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37
	for 22" drive wheel						
	<input type="checkbox"/> 33.5	<input type="checkbox"/> 34.5	<input type="checkbox"/> 35.5	<input type="checkbox"/> 36.5	<input type="checkbox"/> <u>37.5</u>	<input type="checkbox"/> 38.5	<input type="checkbox"/> 39.5
	for 24" drive wheel						
	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> <u>39</u>	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42

Seat height front **SHv**

<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47				
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	--	--	--	--

Tipping point **KP**

(Front of back tube up to
axis middle drive wheel)

<input type="checkbox"/> <u>5</u>	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10						
-----------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------	--	--	--	--	--	--

Please cross on the appropriate dimension (all dimensions are specified in cm).

Underlined values are standard settings.

The values SHh and SHv may deviate a little from the specified values depending on the wheel size, wheel camber as well as the front and rear frame angle!

Frame colour – one-colour design		
Information about current colour designs and example images at: www.proactiv-gmbh.de/farben.php		
frame colour: Selection from colour palette		
Frame colour is specified when ordering <input checked="" type="checkbox"/> A9000100170-000		
Powder coatings		
White RAL 9016 <input checked="" type="checkbox"/> gloss A9000100170-023 <input type="checkbox"/> matt A9000100170-123 € 40.10	Brilliant silver <input checked="" type="checkbox"/> gloss A9000100170-011	Anthracite RAL 7016 <input checked="" type="checkbox"/> gloss A9000100170-001 <input type="checkbox"/> matt A9000100170-101 € 40.10
Black RAL 9005 <input checked="" type="checkbox"/> gloss A9000100170-003 <input type="checkbox"/> matt A9000100170-019 € 40.10	Pure orange RAL 2004 <input checked="" type="checkbox"/> gloss A9000100170-022 <input type="checkbox"/> matt A9000100170-122 € 40.10	Pink magenta RAL 4010 <input checked="" type="checkbox"/> gloss A9000100170-037 <input type="checkbox"/> matt A9000100170-137 € 40.10
Pastel blue RAL 5024 <input checked="" type="checkbox"/> gloss A9000100170-036 <input type="checkbox"/> matt A9000100170-136 € 40.10	Lemon green <input checked="" type="checkbox"/> gloss A9000100170-005 <input type="checkbox"/> matt A9000100170-105 € 40.10	Yellow RAL 1023 <input checked="" type="checkbox"/> gloss A9000100170-046 <input type="checkbox"/> matt A9000100170-146 € 40.10
Red RAL 3001 <input checked="" type="checkbox"/> gloss A9000100170-047 <input type="checkbox"/> matt A9000100170-147 € 40.10	Lagoon blue <input checked="" type="checkbox"/> gloss A9000100170-052 <input type="checkbox"/> matt A9000100170-152 € 40.10	Athletic grey <input checked="" type="checkbox"/> gloss A9000100170-050 <input type="checkbox"/> matt A9000100170-150 € 40.10
Muddy grey <input checked="" type="checkbox"/> gloss A9000100170-051 <input type="checkbox"/> matt A9000100170-151 € 40.10		
“Dormant” powder coatings		
Deep black, silver effect <input checked="" type="checkbox"/> gloss A9000100170-007	Sky blue <input checked="" type="checkbox"/> gloss A9000100170-017	Marine blue <input checked="" type="checkbox"/> gloss A9000100170-045 <input type="checkbox"/> matt A9000100170-145 € 40.10
Electric green <input checked="" type="checkbox"/> gloss A9000100170-015 <input type="checkbox"/> matt A9000100170-115 € 40.10	Violet <input checked="" type="checkbox"/> gloss A9000100170-044	Fire red <input checked="" type="checkbox"/> gloss A9000100170-018 <input type="checkbox"/> matt A9000100170-118 € 40.10
Frame colour – “liquid coating” one-colour design		
€ 323.20		
Storm (ice blue) <input type="checkbox"/> gloss A9000100170-500 <input type="checkbox"/> matt A9000100170-501	Coffee <input type="checkbox"/> gloss A9000100170-520 <input type="checkbox"/> matt A9000100170-521	Prairie (rose gold) <input type="checkbox"/> gloss A9000100170-504 <input type="checkbox"/> matt A9000100170-505
Pearl (mother of pearl-white) <input type="checkbox"/> gloss A9000100170-524 <input type="checkbox"/> matt A9000100170-525	Berry <input type="checkbox"/> gloss A9000100170-526 <input type="checkbox"/> matt A9000100170-527	Blue ocean <input type="checkbox"/> gloss A9000100170-522 <input type="checkbox"/> matt A9000100170-523
Emerald (emerald green) <input type="checkbox"/> gloss A9000100170-514 <input type="checkbox"/> matt A9000100170-515	Mint <input type="checkbox"/> gloss A9000100170-528 <input type="checkbox"/> matt A9000100170-529	Chilli <input type="checkbox"/> gloss A9000100170-516 <input type="checkbox"/> matt A9000100170-517

Frame colour – two colour designs

Note: Frames with two colour designs in the gloss level of the design colour.

- ☐ **Two-colour design comprising one basic colour (on the back of the frame) and one design colour (on the front of the frame)** € 248.30

Please select from the following colours:

Base colour (on the rear frame)

White RAL 9016 <input type="checkbox"/> A9000100170-200	Brilliant silver <input type="checkbox"/> A9000100170-201	Anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
Black RAL 9005 <input type="checkbox"/> A9000100170-203	Pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	Lemon green <input type="checkbox"/> A9000100170-207
Yellow RAL 1023 <input type="checkbox"/> A9000100170-208	Lagoon blue <input type="checkbox"/> A9000100170-212	Athletic grey <input type="checkbox"/> A9000100170-210
Muddy grey <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

Black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	Lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	Chilli <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
Pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	Coffee <input type="checkbox"/> gloss A9000100170-278 <input type="checkbox"/> matt A9000100170-279	Mint <input type="checkbox"/> gloss A9000100170-286 <input type="checkbox"/> matt A9000100170-287
Emerald (emerald green) <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	Storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265	Prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267
Pearl (mother of pearl-white) <input type="checkbox"/> gloss A9000100170-282 <input type="checkbox"/> matt A9000100170-283	Berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	Blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: If no specifications have been made, the caster fork and bearing block are mounted in black.

Drive wheel: Aluminium rim, hub and spokes

- ☐ Black A9000100171-003

Caster fork and bearing block

- ☐ Silver A9000100173-011
☐ Black A9000100173-003

Attachment parts colour design according to customer desires:

Drive wheel: Aluminium rim and design inlet¹

- ☐ Pure orange RAL 2004 A9000100171-022 € 193.10
☐ White RAL 9016 A9000100171-023 € 193.10
☐ Frame colour A9000100171-999 € 230.80
☐ Selection of single colour coating (deviating from the frame colour) A9000100171-998 € 267.10
Colour: _____
A9000100170- ☐☐☐

Caster wheel: 5" aluminium rim² (3-spoke design)

- ☐ Pure orange RAL 2004 A9000100174-022 € 44.40
☐ White RAL 9016 A9000100174-023 € 44.40
☐ Frame colour A9000100174-999 € 69.50
☐ Selection of single colour coating (deviating from the frame colour) A9000100174-998 € 80.60
Colour: _____
A9000100170- ☐☐☐

Special colour requests (frame or attachment parts)

Frame colour individual procurement Selection from RAL colour palette

- ☐ A9000100170-100 € 383.10

Degree of gloss

- ☐ gloss
☐ matt

Colour designation

RAL number

Note: Attachment parts **not** possible in design colours

¹ Only in conjunction with Liteline hub (for drive wheel aluminium rim with 18 or 36 spokes, object).

² Only in conjunction with 5" caster wheel aluminium rim 3-spoke design easy running A9000501200 or good shock absorption A9000501210

Frames

Custom-made seat width 39 cm - 46 cm


☐ A9041100060 € 282.90

with frame reinforcement, double cross-tube bolted, side panel screw connection via the reinforcement inlet

Caster wheel bearing block

Base plate for bearing block welded

☒ A9041500150



Base plate for bearing block clamped

☐ A9041500600-K* € 108.60



Push handles

Safety push handles offset to the rear

☒ Length: 45 cm
A9041401075-K001
Max. height = SHh + approx. 53 cm

☒ Length: 55 cm
A9041401075-K002
max. height = SHh + approx. 63 cm
continuous height adjustment
Safety push handle holder mounted on the back cross bar



Safety push handles aluminium offset to the rear with side handle offset 145 mm

☐ A9041401080-K* € 82.80

continuous height adjustment
Safety push handle holder mounted on the back cross bar
Length: 55 cm
Up to seat width of 38 cm or body weight of 60 kg



Central push handle, aluminium Foldable handle area

☐ A9041401086-K € 94.90


Height-adjustable and foldable in the handle area via toothed joint
Holder welded on the middle on the back cross bar
Up to seat width of 38 cm or body weight of 50 kg
Not compatible with headrest fixture
Note: Not suitable for carrying the wheelchair - to do this, grasp the back cross bar.



Central push handle for brake handle folding gripping area

☐ A9041401070-K € 131.20


Closed, with continuous height adjustment, with foam rubber handle cover.
Up to seat width of 38 cm or body weight of 50 kg
Only in conjunction with parking knee lever brake for accompanying person or drum brake
Not compatible with headrest fixture



Push handle bar aluminium, closed, foldable via toothed joint

☐ A9041401069-K € 360.50

Continuous height adjustable push handle holder back-positioned
Up to seat width of 38 cm or body weight of 60 kg



Backrest

Backrest fixed

☒ A9041400020
bolted with aluminium side sections
Back angle 90° to the ground
cannot be folded



Folding backrest

☐ A9041400010-K € 283.20
7-way inclination adjustment
With integrated aluminium side sections



Back equipment

Adaptable back with back padding

☒ A9041802010-K
Cover: Black, breathable, open surface
Equipment with back cross bar



Ergonomic aluminium back shell from RH 24 cm

☐ A9000802210-K € 408.30
With lumbar support
Colour: Anthracite matt
(Where not specified under colour selection)
Back cushions, if desired, please select separately



Back cushion for ergonomic back shell

☐ A9000802220-K € 201.80
Cushion thickness: 2 cm
Cover:
☐ Black, breathable, open surface
☐ Black, breathable, closed surface
If not specified, the cover fabric "black, breathable, open surface" is delivered.



ERGO Back backrest bracket for ergonomic back shell from RH 24 cm and SB 30 cm

☐ A9000802019-K € 224.30
With recessed gripping area
Colour: Anthracite matt
(Where not specified under colour selection)
Only compatible with central push handle art. no. A9041401082
Back shell and back cushions, if desired, please select separately



Equipment without back upholstery and back padding or back system


☐ A9000802011 € -53.00
Equipment with back cross bar



Seat upholstery

Body Contour

☒ A9000801025-K



Aluminium seat plate

☐ A9000801210-K € 270.70
Assembly between the frame tubes,
on bolted welded seat plate holders
Colour: Anthracite matt
(Where not specified under colour selection)



Equipment without seat upholstery

☒ A9000801050



Colours of back equipment and seat

Information about current colour designs and availability at www.proactiv-gmbh.com

Anthracite RAL 7016 gloss

- ☒ A9000100179-995
☐ Back shell
☐ ERGO Back backrest bar
☐ Aluminium seat plate

in frame colour

- ☐ A9000100179-999 € 85.30
☐ Back shell
☐ ERGO Back backrest bar
☐ Aluminium seat plate

For frame coating in two-colour design

Please specify whether the desired option is to be coated in the basic colour or the design colour

- ☐ A9000100179-999 € 85.30
☐ Back shell
☐ Base colour
☐ Design colour
☐ ERGO Back backrest bar
☐ Base colour
☐ Design colour
☐ Aluminium seat plate
☐ Base colour
☐ Design colour

If another colour is desired, please select this from the wheelchair colour scheme

In colour from one-colour design:



- ☐ A9000100179-998 € 138.10
 Please enter the last 3 digits of the desired colour:
 A9000100170-
☐ Back shell
☐ ERGO Back backrest bar
☐ Aluminium seat plate

In colour from two-colour design:

- ☐ A9000100179-997 € 223.60
 Please enter the last 3 digits of the desired colour:
 A9000100170-
☐ Back shell
☐ ERGO Back backrest bar
☐ Aluminium seat plate

Unless otherwise specified, the seat and back system is supplied in the standard colour of anthracite matt.

Seat cushion



Seat cushion in rectangular contour Inlet: Standard

- ☐ A9000805180-K* € 106.60

Seat cushion with radius contour Inlet: Standard

- ☐ A9000800300-K* € 148.00

Seat cushion height:

- ☐ 3 cm ☐ 5 cm
☐ 4 cm ☐ 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

- ☐ up to 30x30 cm
 A9000805204-K* € 276.20
☐ >30x30 up to
 50x50 cm
 A9000805200-K* € 401.90

Seat cushion with radius contour Inlet: Anti-Decubitus

- ☐ up to 30x30 cm
 A9000800311-K* € 326.50
☐ >30x30 up to
 50x50 cm
 A9000800310-K* € 443.70

Seat cushion height:

- ☐ 4 cm
☐ 6 cm

Cover:

- ☐ Black, breathable, open surface
☐ Black, breathable, closed surface

If not specified, the cover fabric "black, breathable, open surface" is delivered.

☐ Seat cushion is existing

Brand/type _____

Cushion height _____ cm

You can also find the seat cushions on our separate "PRO ACTIV Seat and Back Systems" order sheet

Clothing guard

Clothing guard

- ☐ A9041403110-K € 345.00
 plug-in, removable



Aluminium clothing guard

Permanently mounted on the side guard
Only in combination with fixed backrest

In black gloss

☐ A9000403100-K € 224.90

in frame colour (basic colour with two-colour design)

☐ A9041403101-K € 287.60

in design colour (only possible with frame colour in two-colour design)


☐ A9041403102-K € 319.00

Armrests

Flip-down armrest with arm padding 25 cm

☐ A9041403140-K € 283.40

Not in combination with clothing guard



Arm padding 25 cm length mounted on clothing guard

☐ A9000403151 € 153.30

Only in combination with clothing guard
Depending on the clothing guard adjustment, the upper edge of the arm padding is mounted approx. 20 to 24 cm above the seat tube.



Axle plate

Standard equipment with 12 mm sockets/axle. Further axle plate on request.

2-slot axle plate

☐ A9000201000*
Quick-release axle Ø 12 mm
Rear SH:
24" wheel: approx. 36 - 42 cm
22" wheel: approx. 33.5 - 39.5 cm
20" wheel: approx. 31 - 37 cm




1-slot axle plate

☐ A9000201001*
Only available in black
(no colour selection)
Only in combination with 10° camber
Rear SH can only be adjusted in some cases



Extended bike plate

☐ A9000201030* € 291.40
+ 4° camber integrated
Adjustable track



Wheel camber

*Total width with tight-mounted handrim

1° drive wheel camber

☐ A9000202010*
Total width = approx. SB + 20 cm*

4° drive wheel camber

☐ A9000202030*
Total width = approx. SB + 26 cm*

6° drive wheel camber

☐ A9000202040*
Total width = approx. SB + 30 cm*

8° drive wheel camber

☐ A9000202080* € 71.20
Total width = approx. SB + 33 cm*
Wheel plate mounted on the outside

10° drive wheel camber

☐ A9000202050* € 126.60
Total width = approx. SB + 39 cm*
Wheel plate mounted on the outside

E-adaptations

Preparation for "e-motion" mounting

☐ A9000201200* € 465.00
Equipment: e-motion torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber.

Preparation for "e-fix" mounting

☐ A9000201210* € 465.00
Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber, as well as BUDDY mount.

Wheelchair delivery to

☐ Alber GmbH
(Please also enter on p.1 under "Shipping address, other")

☐ Specialist dealer

Brake system

Brake system push to lock Standard brake handle

☒ A9000701010-001



Brake system push to lock Long brake lever

☒ A9000701010-002



Brake lever at 20°

☒ A9000700002

Only in conjunction with knee lever brake system, standard or long brake lever



Brake system push to lock Brake lever with M8 adapter and plastic ball

☐ A9000706001 **€ 154.70**

Only in conjunction with knee lever brake system, standard or long brake lever

Brake lever:

☐ Standard

☐ Long



Brake system push to lock Brake handle to flip-down

☐ A9000706008 **€ 73.40**

Length: 9 cm



Knee lever brake with one-hand operation

☐ A9000707020-K **€ 362.80**

via Bowden cable

Operation:

☐ Right

☐ Left



Parking knee lever brake for accompanying person

☐ A9000701020 **€ 429.10**

Push handle operation via Bowden cable.



Drum brake

☐ A9000702011-K* **€ 438.30**

Operating from the push handle

Total width + 3cm (see wheel camber)

Please note:

- Only in combination with drive wheel aluminium 36-spoke rim wheel
- Not in combination with one-hand drive, e-fix or e-motion



Knee lever brake with lock

☐ Standard brake handle A9000700214 **€ 146.70**

☐ Long brake lever A9000700215 **€ 146.70**

Locking block for brake lever

Prevents braking during pushing by accompanying persons

Not possible in combination with knee lever brake with one-hand operation



Footrest

Continuous footrest

☒ A9041601010-K*
with tubular resting bracket
height-adjustable
from USL 26 cm to USL 47 cm possible
For equipment of the footrest with "Rubber
coated footplate support" please select under
"Footrest accessories".



Footrest folding up to one side

☐ A9041601300-K* € 160.80
from USL 25 cm:
with continuous tubular resting bracket,
height adjustable +/- 3 cm,
cannot be removed.
from USL 29 and frame height 21 cm:
Can also be pivoted outwards.
For equipment of the footrest with "Rubber
coated footplate support" please select under
"Footrest accessories".



Footbar folds up to the rear with spring locking mechanism

☐ A9041602005-K* € 265.70
with continuous tubular resting bracket
from USL 32 cm
For equipment of the footrest with "Rubber
coated footplate support" please select under
"Footrest accessories".



Divided footrest

☐ A9041603200-K € 311.40
with rubber coated footplate, height-
adjustable, can also be pivoted outwards and
folded up to the side, cannot be removed
from USL 29 cm and SB 26 cm possible



Footrest, positioned on the upper frame

☐ A9041606010-K
height-adjustable
From USL 20 cm to USL 25 cm possible.
Rubber coated footplate support required,
please select this under "Footrest
accessories".



Footrest accessories

Rubber coated footplate, 14 cm deep

☐ A9000610000-K € 116.40
bolted, angle-adjustable



Rubber coated footplate, 18 cm deep

☐ A9000610005-K € 129.50
bolted, angle-adjustable



Footplate side stop

☐ A9000610100 € 55.20
adjustable
Only in combination with rubber coated
footplate.



Deflector rollers mounted on the footplate

☐ A9000600006 € 108.60
Only in combination with rubber coated
footplate.



Adaptations

Anti-tipping support

- ☐ right side
A9041301160-K001 € 184.00
- ☐ left side
A9041301160-K002 € 184.00
- with tension spring, foldable
removable using Quick-pin



Outdoor push bar

- ☐ A9041401210 € 419.50
- for assistance by an accompanying person,
adaptable to the frame cross bar from behind
removable via quick release clips
Possible from SB 19 cm



Outdoor wheel attachment with steering fork

- ☐ A9041900100 € 782.00
- Equipped with straight-running lock
can be removed from the wheelchair via quick
release lever
can be removed from the support tube via
quick-release axle with swivelling pneumatic
tyre 12.1/2 x 2.1/4" caster wheel
for maximum ground pressure at the steering
attachment of 30 kg



Transport restraint connection integrated on the rear frame with belt eye

- ☐ A9040900001 € 337.90
- Frame-welded part for mounting a lap belt and
the rear retractor. Lap belt and buckle
(A9000900001 / A9000900007), please select
separately if required.



Bearing block intermediate plate with belt eye for mounting the front retractor

- ☐ Colour: black
A9000900008-001 € 168.10
- ☐ Colour: Silver
A9000900008-002 € 168.10

Please note:
The wheel track guide increases by
12 mm per side.



Headrest with segment padding

- ☐ A9000900004 € 384.30
- 3-part
and support tube ø 25 mm



Headrest fixture

- ☐ pipe dimensions support
tube ø 25 mm € 258.60
A9000900005-K001
- For mounting the headrest segment padding
Please select separately (A9000900004).
Reinforced version with double cross-tube.
- ☐ pipe dimensions support
tube 15 x 15 mm € 207.60
A9000400850-001
- ☐ pipe dimensions support
tube
12.7 x 12.7 mm € 207.60
A9000400850-002
- with square holding fixture for headrests.
Assembly on the back cross bar.
Reinforced version with double cross-tube.
Possible from SB 23 cm.
Not in combination with the central push
handle.



You will find our PRO ACTIV transport restraint systems
on our separate "Transport restraint systems" order sheet



3" caster wheel

3" Skate wheel (80 mm)

- ☐ A9000501110 € 65.40
puncture-proof



4" caster wheel

4" caster wheel "Standard" solid rubber black

- ☒ A9000501011-004
(82A), puncture-proof



4" Skate wheel (100 mm)

- ☐ A9000501120 € 65.40
puncture-proof



4" safety light-up wheel (100 mm)

- ☐ A9000501220 € 87.60
puncture-proof



4" Frog Legs caster wheel (100 mm)

- ☐ A9000501230 € 196.50
with replaceable rubber 6 mm, puncture-proof
Colour: Black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel "Standard" PU black

- ☒ A9000501036-005
(82A), puncture-proof



5" Caster wheel "Standard" PU grey

- ☒ A9000501031-005
(82A), puncture-proof



5" Caster wheel "2 components"

- ☐ A9000501030 € 52.00
puncture-proof



5" safety light-up wheel (125 mm)

- ☐ A9000501221 € 87.60
puncture-proof



5" Aluminium rim 3-spoke design "easy running"

- ☐ Silver, anodised A9000501200 € 160.50
☐ black A9000501201 € 160.50
With exchangeable rubber ring 3 mm hard, puncture-proof
Incl. hair guard spacer (3 mm)



5" aluminium rim 3-spoke design "Good shock absorption"

- ☐ Silver, anodised A9000501210 € 166.30
☐ black A9000501211 € 166.30
with replaceable rubber 6 mm medium hard, puncture-proof
Incl. hair guard spacer (3 mm)



5" Frog Legs caster wheel (125 mm)

☐ A9000501232 € 196.50

with replaceable rubber 6 mm, puncture-proof
Colour: Black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" Caster wheel "Standard" PU grey

☐ A9000501010-006
(82A), puncture-proof



6" caster wheel "pneumatic tyre"

☐ A9000501050 € 52.90
SH front min. 42 cm



Caster fork

Caster fork colours, see "Component colours"

Caster fork Standard

☐ A9000500100-K001
with screwed axle

☐ A9000500100-K002 € 154.80
with quick-release axle



Heavy duty caster fork

For high loads or
when using additional electric drives

☐ A9000510111-K001 € 49.90
with screwed axle

☐ A9000510111-K002 € 202.40
with quick-release axle



Suspension fork

incl. 5" aluminium rim "good shock absorption"

☐ A9000500160-K001 € 466.90
with screwed axle

☐ A9000500160-K002 € 570.40
with quick-release axle



Caster fork spare wheel Standard

☐ A9000500300-K € 274.40
Standard caster forks with quick-release axles



Caster fork spare wheel Heavy Duty

☐ A9000500301-K € 312.50
Heavy duty caster forks with quick-release
axles



The spare caster fork does not have a caster
wheel.

In case a caster wheel is requested, please
put a cross on a second selection in the
category "Caster wheels".

Track distance

Caster wheel track extension



☐ A9000500400* € 54.70
10 mm per side

☐ A9000500401* € 54.70
15 mm per side

☐ A9000500410* € 54.70
20 mm per side

Drive wheels

Drive wheel, 18 spoke aluminium rim

☒ 20" (ETRTO 451-19)
A9000203052

☒ 22" (ETRTO 489-19)
A9000203062

☒ 24" (ETRTO 540-19)
A9000203122

with Liteline hub, spoked radially

Not in combination with:

- Puncture-proof tyres
- One-hand drive
- Drum brake
- Adaptive bike

only up to a body weight of 70 kg:

Sturz \ Gewicht	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel, 36 spoke aluminium rim

☒ 20" (ETRTO 451-19)
A9000203050

☒ 22" (ETRTO 489-19)
A9000203060

☒ 24" (ETRTO 540-19)
A9000203000

Up to 120 kg body weight with Liteline hub (left image)

Over 120 kg body weight with aluminium high flange (right image)
radially spoked




Drive wheel double rim 36 spokes

☐ 24" (ETRTO 540)
A9000203100

€ 145.90

with aluminium front hub,
radially spoked
Use with high loads



Spinergy running wheel "Spox Sport" with 24 PBO spokes Rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203550

€ 725.00



Spinergy "Everyday Light Extreme" with 12 PBO spokes Rim and hub in black, spokes in white

☐ 24" (ETRTO 540)
A9000203554

€ 885.00



Spinergy "Everyday Light Extreme" with 12 PBO spokes Rim, hub and spokes in black


☐ 24" (ETRTO 540)
A9000203552

€ 885.00



Handrims

x = dimension handrim attachment



Handrim attachment

General:

☐ x = 11 mm

☐ x = 19 mm

Handrim attachment

Curve handrim:

☐ x = 10 mm

☐ x = 18 mm

Unless otherwise specified, the handrim is mounted closely (11 mm / 10 mm)

Stainless steel Ø 20 mm

☐ 20" handrim for
20" wheel
A9000205010-001

☐ 22" handrim for
22" wheel
A9000205010-002

☐ 24" handrim for
24" wheel
A9000205010-004



Aluminium Ø 19 mm silver anodised

☐ 24" handrim for
24" wheel
A9000205019-004



Aluminium Ø 16 mm powder coated in the frame colour

(basic colour with two-colour design)

☐ = Gloss surface

☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

☐ 20" handrim for
20" wheel
A9000205021-001 **€ 63.10**

☐ 22" handrim for
22" wheel
A9000205021-002 **€ 63.10**



Aluminium ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

☐ = Gloss surface

☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

☐ 20" handrim for
20" wheel
A9000205020-001 **€ 63.10**

☐ 22" handrim for
22" wheel
A9000205020-002 **€ 63.10**

☐ 24" handrim for
24" wheel
A9000205020-004 **€ 63.10**



Aluminium ø 16 mm polished and hard-anodised

☐ 20" handrim for
20" wheel
A9000205031-001 **€ 152.50**

☐ 22" handrim for
22" wheel
A9000205031-002 **€ 152.50**




Aluminium ø 19 mm polished and hard-anodised

☐ 20" handrim for
20" wheel
A9000205030-001 **€ 152.50**

☐ 22" handrim for
22" wheel
A9000205030-002 **€ 152.50**

☐ 24" handrim for
24" wheel
A9000205030-004 **€ 152.50**



Titanium ø 19 mm <input type="checkbox"/> 24" handrim for 24" wheel A9000205040-004 € 310.30 	Silicone handrim cover <input type="checkbox"/> for 20" handrim A9000205100-001 € 215.20 <input type="checkbox"/> for 22" handrim A9000205100-002 € 215.20 <input type="checkbox"/> for 24" handrim A9000205100-004 € 215.20 Surface without structure compatible with handrim Ø 19 mm or Ø 20 mm 	Ultra-Grip "very grippy" <input type="checkbox"/> 20" handrim for 20" wheel A9000205070-001 € 372.10 <input type="checkbox"/> 22" handrim for 22" wheel A9000205070-002 € 372.10 <input type="checkbox"/> 24" handrim for 24" wheel A9000205070-004 € 372.10 
Ergo-Grip <input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 377.90 <input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 377.90 <input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 377.90 Attachment fitting to the rim 	Maxgrepp with cell and natural rubber cover very good grip, soft <input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 377.90 <input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 377.90 <input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 377.90 	Curve black handrim anodised, ergonomic handrim profile <input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205130-004 € 370.60 <input type="checkbox"/> 24" handrim for 24" wheel H 31 mm size L A9000205140-004 € 370.60 
Curve silver anodised ergonomic handrim profile <input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel H 31 mm, size L A9000205160-004 € 375.50 	Curve powder-coated, ergonomic handrim profile <input type="checkbox"/> = Frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface) <input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205135-004 € 432.40 <input type="checkbox"/> 24" handrim for 24" wheel H 31 mm size L A9000205145-004 € 432.40 	Curve Grip, with Tetra-Grip coating ergonomic handrim profile <input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205155-004 € 481.40 

<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004 € 492.90</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205175-004 € 477.30</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205185-004 € 488.80</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
<p>Gekko S black anodised</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205195-002 € 477.30</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205195-004 € 477.30</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	
<p>Spoke guard</p> <p>Finger- and spoke guard "transparent"</p> <p><input type="checkbox"/> for 20" wheel A9000206300-001 € 191.50</p> <p><input type="checkbox"/> for 22" wheel A9000206300-002 € 191.50</p> <p><input type="checkbox"/> for 24" wheel A9000206300-004 € 191.50</p> 		

Tyres

Schwalbe "RightRun" High-pressure tyre

- ☐ 22x1" (489-25)
A9000204001-002
- ☐ 24x1" (540-25)
A9000204001-004
- Colour: black
7.5 bar



Schwalbe "RightRun" High-pressure tyre

- ☐ 20x1" (451-25)
A9000204000-001
- ☐ 22x1" (489-25)
A9000204000-002
- ☐ 24x1" (540-25)
A9000204000-004
- Colour: Grey, black wall
7.5 bar



Schwalbe "DownTown" Profiled tyres

- ☐ 24x1" (540-25)
A9000204010-004
- 7.5 bar



Schwalbe "MarathonPlus Evolution" with "smart guard" puncture protection and reflector strips

- ☐ 22x1" (489-25) € 45.50
A9000204030-002
- ☐ 24x1" (540-25) € 45.50
A9000204031-004
- Colour: black
10 bar



Schwalbe "MarathonPlus Evolution" with "smart guard" puncture protection

- ☐ 24x1" (540-25) € 45.50
A9000204030-014
- Colour: black
10 bar



Schwalbe "MarathonPlus" with "smart guard" puncture protection

- ☐ 24x1" (540-25) € 45.50
A9000204030-004
- Colour: Grey, black wall
10 bar



Schwalbe "One"

- ☐ 24x1" (540-25) € 86.00
A9000204123-004
- Colour: black
9 bar



Puncture-proof tyre "solid rubber tyres"

- ☐ 24x1" (540-25) € 113.40
A9000204201-004
- Colour: black
Not in combination with the 18 spoke drive
wheel, Spinergy and one-hand drive



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminium front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

- ☐ A9000203410 € 650.40
24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel













- ☐ A9000203413 € 650.40
24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel



NOTE: We charge a surcharge for the
use of drum brakes

- ☐ A9000200004 € 210.60

Through axle with assistance, Large push button <input type="checkbox"/> A9000200096 € 97.80 Interlock dimension: 122 / 110 mm Axle diameter: 12 mm Not in conjunction with lockable drive wheel axle clip "Tetra Clip", Not in combination with one-hand drive 	Lockable drive wheel axle clip "Tetra Clip" <input type="checkbox"/> A9000206400 € 153.90 Not in combination with Spinerger drive wheels, one-hand drive, drum brake or Liteline hub 	Hand air compressor with pressure gauge <input type="checkbox"/> A9000901100 € 98.30 up to 14 bar 
Length-adjustable safety belt with buckle <input type="checkbox"/> A9000902100 € 74.60 	Adjustable safety belt with hook and fastener <input type="checkbox"/> A9000902110 € 59.60 	Crutch holder with nylon cup <input type="checkbox"/> A9000902310 € 127.00 flexible 
SIGMA lighting set <input type="checkbox"/> A9000903270 € 99.80 Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations 	Transport bag for drive wheels <input type="checkbox"/> for 24" wheel Ø 76 cm x 2 cm A9000903530-004 € 71.90 	Accessory bag <input type="checkbox"/> A9000903600 € 70.30 17 cm x 13cm x 5 cm 
Therapy table <input type="checkbox"/> Width 350 mm (for SB up to 37 cm) € 444.00 A9000902800-001 <input type="checkbox"/> Width 400 mm (SB 38 - 41 cm) € 444.00 A9000902800-002 <input type="checkbox"/> Width 420 mm (from SB 42 cm) € 444.00 A9000902800-003 removable table plate from Lexan, transparent, with front edging, with 2 items Positioning stops Ø 25 mm (for adaptation to frame tube diameter Ø 25 mm) 		
Further accessories, e.g., adaptive bikes, steering attachments, etc. see separate order sheet or www.proactiv-gmbh.com		

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 20 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs ☐ = Included in the basic price ☒ = Position selected

Order sheet BUDDY classic Rev. 07 UK