PRO ACTIV Reha-Technik GmbH
Im Hofstätt 11
D-72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-mail: info@proactiv-gmbh.de
www.proactiv-gmbh.com



PED-E € 4,142.02

- Aluminium frame with CNC milled/welded construction
- Aluminium fork with milled fork bridges,
  Colour: black
- Hydraulic disc brakes 180 mm / 160 mm
- Neodrives Z15 drive system with TFT display
- Rechargeable battery integrated in the frame
  37 V; 16 Ah; 587 Wh
- Shimano Zee dérailleur 1x10
- Pedal crank (L = 80 mm)
- Aluminium rims
- Saddle support Ø 27.2 mm
- Side stand
- Bell
- Lighting set
- Cargo load: 75 kg

☐ **Offer date:** _____

☐ **Order date:** _____
according to GTC PRO ACTIV Reha-Technik GmbH

Consultant: _____

Order address / customer:

Cust. no.: _____

Last name: _____

First name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-mail: _____

Collection: ☐ Yes ☐ No

Delivery address

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Date of birth: _____

Note: _____

Biometric data:

Body height: _____

Weight: _____

Inside leg: _____

◧ = Option included in basic price ☐ = Option with surcharge; please select the desired options: ☒ ; ☒

Wheel size

18"

◧ A9525000180

18x2" Schwalbe Big Apple

◧ A9525200010
(50-355)
Performance Line
RaceGuard E25
Puncture protection level 4
Colour: black



20"

◧ A9525000200

20x2.35" Schwalbe "Smart Sam"

☐ A9525105101 **€ 24.79**
(60-406)
Performance Line
Puncture protection level 3.5
Colour: black



20x2.15" Schwalbe "Big Ben Plus"

☐ A9525105102 **€ 24.79**
(55-406)
Performance Line
DD, GreenGuard E50
Puncture protection level 6
Colour: black



Bottom bracket position

Please specify the desired bottom bracket position:

_____

If nothing is specified, the standard bottom bracket position is 1.



Saddle

Velo Fit Traveller saddle

◧ A9500800010
Size M
Width: 15 cm; length: 27 cm

◧ A9500800020
Size L
Width: 19 cm; length: 26 cm

Colour: black



Definition of the saddle support length

Distance from the upper edge of the rechargeable battery to the upper edge of the saddle:

_____ cm

If nothing is specified, the saddle support will be delivered with extended length (40 cm) and can be adjusted by the customers themselves.



◧ = Option included in basic price ☐ = Option with surcharge; please select the desired options: ☒ ; ☒

Front end

Front end 7°

◧ Length: 90 mm
A9500400500-090

◧ Length: 110 mm
A9500400500-110

Clamping: Oversize 31.8 mm
Colour: black



Ergotec Swell RE front end

☐ A9500400550 **€ 39.58**

Length: 100 mm
Clamping: Oversize 31.8 mm
angle-adjustable
Colour: black



Handlebar

Sport Ergotec handlebar

◧ A9500400610
Width: 720 mm
Height: 15 mm
Grip bend: 9°
Oversize 31.8 mm
Colour: black

The handlebar width can be individually adjusted in the range from 720 to 615 mm: _____ mm



Comfort Ergotec handlebar

◧ A9500400620
Width: 680 mm
Height: 30 mm
Grip bend: 14°
Oversize 31.8 mm
Colour: black

The handlebar width can be individually adjusted in the range from 680 to 615 mm: _____ mm



Grips

Sport grips black / grey

◧ A9500400700
Shock-absorbing gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing grip cover

☐ A9500400710 **€ 21.26**

Ergonomic gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing Comfort grip cover

☐ A9500400720 **€ 26.81**

Large, ergonomically shaped wrist rest,
anti-twist clamping on the handlebar



Pedals

Offroad pedals

◧ A9500500300
CR-MO axle,
with reflector (K-mark)
Platform with jagged profile



Touring pedals

◧ A9500500310
CR-MO axle,
with reflector (K-mark),
rubberised support



Ergo Grip pedals

☐ A9500500320 **€ 26.13**

CR-MO axle,
with reflector (K-mark),
sandpaper support



◧ = Option included in basic price ☐ = Option with surcharge; please select the desired options: ☒ ; ☒

Accessories

Mudguard 53 mm,
Colour-coated on front and back

☐ A9525700053 **€ 191.43**

not for 16"
Colour selection:
see colour design at the back




Luggage rack

☐ A9525700100 **€ 416.47**

Aluminium
Incl. rear light with battery



Colour design – frame and attachment parts

For the following components, a colour can be chosen from the information sheet "Pedelec Colour Design":

RAL powder paint or "Dormant" powder paint or "Liquid coating" custom paint

Frame colour:

◧ RAL powder paint and "Dormant" powder paint:

Please enter the last 3 digits of the colour:

A9500100170 - ☐☐☐ in words (colour): _____

☐ "Liquid coating" custom paint: **€ 161.68**

Please enter the last 2 digits of the custom colour:

A9500100170 - 5☐☐ in words (colour): _____

Mudguard colour:

◧ RAL powder paint or "Dormant" powder paint:

A9500700170

Please enter the last 3 digits of the colour:

A9500100170 - ☐☐☐ in words (colour): _____

☐ "Liquid coating" custom paint:

A9500700172 **€ 111.85**

Please enter the last 2 digits of the custom colour:

A9500100170 - 5☐☐ in words (colour): _____

Notes and remarks:

I hereby confirm the information provided above and agree to receiving a written offer without obligation.

Location, date	**End customer signature (offer)**	**Location, date**	**Consultant signature (offer)**
Location, date	**End customer signature (order)**	**Location, date**	**Consultant signature (order)**

◧ = Option included in basic price ☐ = Option with surcharge; please select the desired options: ☒ ; ☒

Info sheet: Colour Designs Pedelecs

Colour Designs - Powder Coatings

Information on current colour designs and sample images at: www.proactiv-gmbh.com/colours/

Powder Coatings

white RAL 9016 ◧ gloss A9500100170-023	**anthracite RAL 7016** ◧ gloss A9500100170-001 ◧ matt A9500100170-101	**black RAL 9005** ◧ gloss A9500100170-003 ◧ matt A9500100170-019
pure orange RAL 2004 ◧ gloss A9500100170-022 ◧ matt A9500100170-122	**yellow RAL 1023** ◧ gloss A9500100170-046 ◧ matt A9500100170-146	**red RAL 3001** ◧ gloss A9500100170-047 ◧ matt A9500100170-147

„Dormant“ Powder Coatings

sky blue ◧ gloss A9500100170-017	**electric green** ◧ gloss A9500100170-015	**violet** ◧ gloss A9500100170-044
fire red ◧ gloss A9500100170-018		

Powder Coatings “Trend colours”

athletic grey ☐ gloss A9500100170-303 **€ 30.59** ☐ matt A9500100170-304 **€ 30.59**	**muddy grey** ☐ gloss A9500100170-305 **€ 30.59** ☐ matt A9500100170-306 **€ 30.59**	**lagoon blue** ☐ gloss A9500100170-301 **€ 30.59** ☐ matt A9500100170-302 **€ 30.59**
peach RAL 3022 ☐ gloss A9500100170-038 **€ 30.59** ☐ matt A9500100170-138 **€ 30.59**		

Colour Designs - Custom Colours ´Liquid Coating´

€ 161.68

storm (ice blue) ☐ gloss A9500100170-500 ☐ matt A9500100170-501	**prairie (rose gold)** ☐ gloss A9500100170-504 ☐ matt A9500100170-505	**blue ocean** ☐ gloss A9500100170-522 ☐ matt A9500100170-523

◧ = Option included in the basic price ☐ = Option with surcharge; please select the desired options: ◧; ☒